

CellSonic VET

VIPP Therapy System

NON-INVASIVE | NO DRUGS | NO SIDE EFFECTS | CHEAPEST COST OF CURE

We tried it on humans first!

A healing machine for animals bears a lot of resemblance to machines used on humans with important differences. An animal cannot tell us what they feel. We have to be careful. Where exactly has the machine to be aimed? What is hurting the animal?


Our experience goes long back when the idea of a machine doing a better job than drugs and a scalpel was ridiculed. The acceptance since then comes from seeing no side effects, a great cost saving and better healing.

A machine for vets has to be rugged. Animals may be unpredictable. Attention is on the animal, not the machine so its controls must be instinctive. If the operator jumps back, the machine must stop on its own.

CellSonic Vet is built for the stables, the yard, an environment far removed from the hygiene of a hospital and yet it delivers the same tissue improving pulses. Bone and skin is made to grow. Nerves re-form. Infection is killed without anti-biotics. Lameness in a horse can be corrected. Injuries on a dog can be healed without disturbing their digestion with drugs.

CellSonic's clever secret is the rapid pulse. Over the years this has been made faster so that it kills germs better and triggers the formation of osteoblasts and fibroblasts. There are easier engineering methods to save production costs but they do not deliver healing so they are not used by CellSonic. By any calculation, CellSonic gives you the lowest cost of cure.

If you have an animal problem not on the standard list, send us your question. What you know for sure is that the treatment is safe- we tried it on humans first.



The CellSonic VIPP machine is a new branch of veterinary medicine that kills infection and triggers tissue regeneration.

The origins of the technology go back forty years to removing human kidney stones by shattering them with sound waves and this was the first ever non-invasive surgery. Since then the technology has been developed with new applications discovered.

Drugs are not used so there are no side effects. VIPP is Very Intense Pressure Pulses which are fast bangs made by flashing high voltages across an electrode in a reflector containing water in a shock head. Pressure wave travels through the body without harming healthy cells. With this Impossible wounds are healed. The infection is killed and the immune system triggered to cause new skin to grow.

CellSonic VET

VIPP Therapy System

NON-INVASIVE | NO DRUGS | NO SIDE EFFECTS | CHEAPEST COST OF CURE



From the healing of broken bones and non-union fractures came the physiotherapy applications and from that the ability to heal and remove the cause of pain. Treatment of muscles increases their strength and this applies to animals involved in sporting activities as well as the infirm. Interestingly, race horses have been cured of lameness and given the strength to win races.

In every case, it can be shown that CellSonic gives the lowest cost of cure and that is the only cost that matters. Many clinical papers are available. It is not CellSonic's policy to pay for papers to be written; they are all from customers and unbiased.

Portability

The complete device fits in a car boot. It can be ported anywhere easily.



Manufactured by :



Apex Meditech

S No.635/1, Plot 49, Vaibhav Society,
Canara Bank Lane, Off Swami Vivekanand Road,
Bibwewadi, Pune 411037. India
Tel : ++91(20) - 24410630, 24410634
E-mail : sales@apex-meditech.com

Quality approved to :
ISO 13485:2016
ISO 9001:2015



www.Apex-Meditech.com

Specifications are subject to change without prior notice