

Micturition Accessories

For Uroflometry



Micturition Chair

The Micturition Chair can be used with the digiFLOW series of UROFLOMETRY Systems. It is a compact and simple chair. It can be used by males as well as females. It is supplied in a dis-assembled form. It can be assembled in a few minutes.

Technical Specifications

Particulars

Dimensions	:	46X38X51 (H) CM
Weight	:	4.6kG
Suitable for	:	Male & Female
Material	:	Frame made with Epoxy powder coated steel, Seat and funnel made of washable plastic
Packing	:	53x46x32 cm., corrugated paper case.

Micturition Accessories

For Uroflometry



Micturition Stand SS

The Micturition Stand can be used with the digiFLOW series of UROFLOMETRY Systems. It is a compact and stand. Good for prostate camps. It can be used by males. It is supplied in a dis-assembled form. It can be assembled in a few minutes. The height of the funnel is adjustable. It is made of Stainless Steel.

Note: The collection jar and the sensor unit are not parts of the Micturition Stand

Technical Specifications

Particulars

Dimensions	:	32X32X53 (H) CM
Weight	:	1.4kG
Suitable for	:	Male
Material	:	Stainless Steel 304
		Funnel made of washable plastic
Packing	:	40x40x35 cm., corrugated paper case.

Micturition Accessories

For Uroflometry

Micturition Stand



The Micturition Stand can be used with the digiFLOW series of UROFLOMETRY Systems. It is a compact and stand. Good for prostate camps. It can be used by males. It is supplied in a dis-assembled form. It can be assembled in a few minutes. The height of the funnel is adjustable.

Technical Specifications

Particulars

Dimensions	:	32X50X64 (H) CM
Weight	:	1.6kG
Suitable for	:	Male
Material	:	Frame made with Epoxy powder coated steel, Funnel made of washable plastic
Packing	:	40x40x35 cm., corrugated paper case.

Quality approved to :
ISO 13485:2016
ISO 9001:2015

Manufactured by :



Apex MediTech

S No.635/1, Plot 49, Vaibhav Society,
Canara Bank Lane, Off S. Vivekanand Road,
Bibwewadi, Pune 411037, India.
Phone. : +91 (20) 24410630, 24410634.
E-mail : sales@apex-meditech.com

www.Apex-Meditech.com