



Optical Fibre Cable

For LED, Xenon and halogen cold light sources
for endoscopies

The illumination guiding you through your
endoscopic procedures guarantees Visual Satisfaction



Optical Fibre Cables for different kinds of Light Sources provide perfect illumination for all endoscopic specialties. Super flexibility of cables makes the cables easy to use. The state-of-the-art fibre material ensures smooth and precise light distribution. Advanced technology optics deliver the bright and balanced illumination physicians and surgeons demand. Excellent engineering, efficient design, and outstanding performance - the qualities you expect from the leader in endoscopic imaging. Available in Ø3.5mm and Ø4.5mm fibre. Standard length is 2.28mtr (7 ½ ft.).

Manufactured by :



Apex MediTech

S No.635/1, Plot 49, Vaibhav Society,
Canara Bank Lane, Off Swami Vivekanand Road,
Bibwewadi, Pune 411037. India
Tel : ++91(20) - 24410630, 24410634
E-mail : sales@apex-meditech.com

Quality approved to :
ISO 13485:2016
ISO 9001:2015



www.Apex-Meditech.com

Specifications are subject to change without prior notice

Optical Fibre Cable