



TIME MACHINE

DIGITAL SKINCARE SYSTEM



SALIENT FEATURES

- Reduces wrinkles considerably, fairly quickly
- Improves skin tone
- Easy to assemble and operate
- 10 Energy Levels
- Digital Energy Level Adjustment
- Time setting from 10 to 20 minutes according to choice
- Four different Phases specially designed to take care of the skin
- Long Lasting Gloves, no need to use sticky or rigid electrodes on face
- Based on an Ultra High Speed Microprocessor
- Two different modes and programmes for Face and Body
- Simple Feather Touch Keys
- Fuzzy Logic for Improved Energy Control
- Portable, Light Weight system
- Free Carrying Case
- Very comfortable treatment
- Unlike conventional TENS no sensation of heavy shocks
- Machine is totally maintenance free
- Consumes low power
- Very low running cost / per treatment cost

TIME MACHINE



TIME MACHINE Technical Specifications

Control Unit

Method	Based on Ultra High Speed Micro Processor
Supply Voltage	91-256V AC @ 50-60 Hz
Power	10 Watts (max)
Fuse	500mA, 5X20mm, Fast Blow
EN60601-1 Class	II
EN60601-1 Type	BF
Power entry	Panel mounting IEC C14 inlet for three wire detachable cable.
Treatment Cycle Timing	10 to 20 minutes programmable with an increment of 1 minute (programmable through Key on front)
Energy Levels	1 to 10, Total Ten choices (programmable through Key on front)
Phases per treatment	Four
Display	Digital Seven Segment LED
Operating temperature	15° - 35° C
Dimensions	265 mm W X 112 mm H X 135 mm D
Weight	2 kG
Degree of protection	IP20

Copper Gloves

Material	Blend of Polymer & Copper
Connection	4 mm button connector
Weight	10 gms (approx.)
Sterilization methods	Ethylene Oxide, Formaldehyde, Gluteraldehyde 2% or Disinfectant Cleaner

Manufactured by :



Apex MediTech

S No.635/1, Plot 49, Vaibhav Society,
Canara Bank Lane, Off Swami Vivekanand Road,
Bibwewadi, Pune 411037. India
Tel : ++91(20) - 24410630, 24410634
E-mail : sales@apex-meditech.com

Quality approved to :
ISO 13485:2016
ISO 9001:2015

www.Apex-Meditech.com

Specifications are subject to change without prior notice