

Trolley



Available in three different tray sizes

Note: The device shown on the Device Tray is not a part of the trolley

www.Apex-Meditech.com

Specifications are subject to change without prior notice



Technical Specifications - Trolley

Type A

Overall Dimensions	870 mm H X 520 mm W X 520 mm D
Tray dimensions	300 mm W X 260 mm D
Weight	10 kG
Gas bottle holder	Ø 100 mm X 200 mm H
Spares compartment	145 mm W X 95 mm D X 75 mm H
Material	Mild steel
Surface finish	Epoxy powder coated

Type B

Overall Dimensions	870 mm H X 520 mm W X 520 mm D
Tray dimensions	270 mm W X 240 mm D
Weight	10 kG
Gas bottle holder	Ø 100 mm X 200 mm H
Spares compartment	145 mm W X 95 mm D X 75 mm H
Material	Mild steel
Surface finish	Epoxy powder coated

Type C

Overall Dimensions	870 mm H X 520 mm W X 520 mm D
Tray dimensions	300 mm W X 180 mm D
Weight	10 kG
Gas bottle holder	Ø 100 mm X 200 mm H
Spares compartment	145 mm W X 95 mm D X 75 mm H
Material	Mild steel
Surface finish	Epoxy powder coated

Manufactured by :



Apex MediTech

S No.635/1, Plot 49, Vaibhav Society,
Canara Bank Lane, Off Swami Vivekanand Road,
Bibwewadi, Pune 411037. India
Tel : ++91(20) - 24410630, 24410634
E-mail : sales@apex-meditech.com

Quality approved to :
ISO 13485:2016
ISO 9001:2015

www.Apex-Meditech.com

Specifications are subject to change without prior notice

Trolley